



NQS6EP

**Electrophoresis
6 fractions**

Description

Sent sample: 1 vial **1 ml liquid sample**.

Storage:

Store the sample at **2-8 °C**.

Instruction for use:

Bring the sample to room temperature at least 1 hour before running the test.

Analyze the control material **as a patient serum**.

The sample is stable **8 hours** at 25 °C, or **until expiration date** if it is kept at **2-8 °C**

Precautions:

The serum has been tested for Human Immunodeficiency Virus (HIV) antibody, Hepatitis B Surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE. However this material should be handled carefully.

Parameters:

Total protein, Albumin, Alfa 1-globuline, Alfa 2-globuline, Beta 1-globuline, Beta 2 -globuline, Gamma-globuline