



NQSSE_P

**Serology
plus**

Description

Sent sample: 1 vial **5 ml lyophilized serum**

Instruction for use:

Bring the sample to room temperature.

Open the vial carefully, and reconstitute the sample with the exact quantity of distilled water (**5 ml**). Leave the sample **30 minute** until the lyophilized material is totally dissolved. During the reconstitution rotate and stir gently the vial for homogenization. Avoid foam formation.

Analyze the control material as a patient sample.

After reconstitution the sample is stable **8 hours** at **25 °C**, or **7 days** if it is kept at **2-8 °C**

Storage:

Store the sample at **2-8 °C**.

Precautions:

The serum has been tested for Human Immunodeficiency Virus (HIV) antibody. However this material should be handled carefully.

Parameters:

HAV IgG, HAV IgM, HBsAg, HCV, HIV, Hiv 1&2, Syphilis (VDRL, RPR, TPHA), HbsAc, Helicobacter Pylori Antibody, Chlamydia IgM and IgG, Cytomegalovirus IgM and IgG, EBV IgM and IgG, Herpes Virus I IgG, Herpes Virus II IgG, Rubella IgG and Rubella IgM, Toxoplasma IgG and IgM, HBcAb, HBcAg, HBeAb, HBeAg