



**NQSPS**

**Specific  
proteins**

**Description**

Sent sample: **1 vial of 5 ml lyophilized sample**

**Instruction for use:**

**1 vial of 5 ml:**

Reconstitute the samples with **5 ml distilled water**.

**Sample reconstitution**

Before reconstitution of the samples, bring it to room temperature.

Open the vial carefully, and reconstitute the sample with the exact quantity of distilled water.

Leave the sample **30 minute** until the lyophilized material is totally dissolved. During the reconstitution rotate and stir gently the vial for homogenization. Avoid foam formation.

Analyze the control material **as a patient serum**.

**Storage:**

Store the sample at **2-8 °C**.

After opening the sample is stable 7 days kept at 25 °C, **2 weeks** kept at **2-8 °C**, **4 weeks** freezed at **-20 °C**

**Precautions:**

The serum has been tested for Human Immunodeficiency Virus (HIV) antibody, Hepatitis B Surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE. However this material should be handled carefully.

**Parameters:**

ASO, CRP, RF, IgA, IgG, IgM,  $\alpha$ 1-antitrypsin,  $\alpha$ 1-glycoprotein,  $\alpha$ 2-macroglobulin, complement C3, complement C4, ceruloplasmin, ferritin, haptoglobin, transferrin