



Instruments:		Rayto RT-7600s, RT-7200, HA-22, Linear STEL 3, SFRI Countender 20 & 30					
Parameter		CONTROL "Low"		CONTROL "Normal"		CONTROL "High"	
		LOT: 1Z0601 L		LOT: 1Z0601 N		LOT: 1Z0601 H	
		Mean ± Limit	Range Variation	Mean ± Limit	Range Variation	Mean ± Limit	Range Variation
WBC	10 ⁹ /L	2,42 ± 0,50	1,92 - 2,92	6,99 ± 1,00	5,99 - 7,99	14,75 ± 2,50	12,25 - 17,25
RBC	10 ¹² /L	2,39 ± 0,20	2,19 - 2,59	4,20 ± 0,30	3,90 - 4,50	5,34 ± 0,40	4,94 - 5,74
HGB	g/dL	7,0 ± 0,5	6,5 - 7,5	12,7 ± 0,6	12,1 - 13,3	16,8 ± 0,8	16,0 - 17,6
HCT	%	21,0 ± 3,0	18,0 - 24,0	37,5 ± 4,0	33,5 - 41,5	47,9 ± 4,5	43,4 - 52,4
MCV	fL	87,9 ± 6,0	81,9 - 93,9	89,2 ± 6,0	83,2 - 95,2	89,7 ± 6,0	83,7 - 95,7
MCH	pg	29,3 ± 3,8	25,5 - 33,1	30,2 ± 3,8	26,4 - 34,0	31,5 ± 3,8	27,7 - 35,3
MCHC	g/dL	33,3 ± 4,8	28,5 - 38,1	33,9 ± 4,8	29,1 - 38,7	35,1 ± 4,8	30,3 - 39,9
RDW-CV	%	14,2 ± 5,0	9,2 - 19,2	14,5 ± 5,0	9,5 - 19,5	17,4 ± 5,0	12,4 - 22,4
PLT	10 ⁹ /L	98 ± 30	68 - 128	241 ± 45	196 - 286	404 ± 65	339 - 469
MPV	fL	7,01 ± 3,00	4,01 - 10,01	7,53 ± 3,00	4,53 - 10,54	7,90 ± 3,00	4,90 - 10,90
PDW	%	14,16 ± 4,00	10,16 - 18,16	17,29 ± 4,00	13,34 - 21,24	18,65 ± 3,00	15,65 - 21,70
PCT	%	0,069 ± 0,050	0,019 - 0,119	0,181 ± 0,100	0,081 - 0,281	0,319 ± 0,200	0,119 - 0,519
LYMPH#	10 ⁹ /L	0,98 ± 0,30	0,68 - 1,28	2,11 ± 0,73	1,38 - 2,84	3,59 ± 1,97	1,62 - 5,56
MID#	10 ⁹ /L	0,17 ± 0,16	0,01 - 0,33	0,50 ± 0,41	0,09 - 0,91	1,12 ± 0,83	0,29 - 1,95
GRAN#	10 ⁹ /L	1,27 ± 0,50	0,77 - 1,77	4,38 ± 0,66	3,72 - 5,04	10,04 ± 3,01	7,03 - 13,05
LYMPH%	%	40,5 ± 10,2	30,3 - 50,7	30,3 ± 10,6	19,7 - 40,9	24,3 ± 8,5	15,8 - 32,8
MID%	%	7,1 ± 5,3	1,8 - 12,4	7,1 ± 5,6	1,5 - 12,7	7,6 ± 6,1	1,5 - 13,7
GRAN%	%	52,4 ± 13,1	39,3 - 65,5	62,6 ± 12,5	50,1 - 75,1	68,1 ± 10,2	57,9 - 78,3

