



Instruments:		DYMIND DH20, DH22, DH26					
Parameter		CONTROL "Low"		CONTROL "Normal"		CONTROL "High"	
		LOT: 1Z0901 L		LOT: 1Z0901 N		LOT: 1Z0901 H	
		Mean ± Limit	Range Variation	Mean ± Limit	Range Variation	Mean ± Limit	Range Variation
WBC	10 ⁹ /L	2,34 ± 0,50	1,84 - 2,84	7,15 ± 1,00	6,15 - 8,15	16,88 ± 2,50	14,38 - 19,38
LYMPH%	%	44,2 ± 8,8	35,4 - 53,0	31,7 ± 7,9	23,8 - 40	29,1 ± 16,0	13,1 - 45,2
GRAN%	%	49,0 ± 17,2	31,9 - 66,2	61,9 ± 9,3	52,6 - 71,2	65,8 ± 8,6	57,3 - 74,4
MID%	%	7,0 ± 5,9	1,0 - 12,9	6,4 ± 5,6	0,8 - 11,9	4,9 ± 4,3	0,6 - 9,2
LYMPH#	10 ⁹ /L	1,04 ± 0,21	0,83 - 1,24	2,27 ± 0,57	1,70 - 2,84	4,92 ± 2,70	2,21 - 7,62
GRAN#	10 ⁹ /L	1,15 ± 0,40	0,75 - 1,55	4,43 ± 0,66	3,76 - 5,09	11,11 ± 1,44	9,67 - 12,56
MID#	10 ⁹ /L	0,16 ± 0,14	0,02 - 0,30	0,45 ± 0,40	0,06 - 0,85	0,83 ± 0,73	0,10 - 1,55
RBC	10 ¹² /L	2,51 ± 0,20	2,31 - 2,71	4,58 ± 0,30	4,28 - 4,88	5,73 ± 0,40	5,33 - 6,13
HGB	g/L	67 ± 6	61 - 73	134 ± 7	127 - 141	172 ± 8	164 - 180
	g/dL	6,7 ± 0,6	6,1 - 7,3	13,4 ± 0,7	12,7 - 14,1	17,2 ± 0,8	16,4 - 18,0
HCT	%	21,8 ± 3,0	18,8 - 24,8	41,9 ± 4,0	37,9 - 45,9	52,5 ± 4,5	48,0 - 57,0
MCV	fL	87,0 ± 6,0	81,0 - 93,0	91,4 ± 6,0	85,4 - 97,4	91,7 ± 6,0	85,7 - 97,7
MCH	pg	26,7 ± 3,8	22,9 - 30,5	29,3 ± 3,8	25,5 - 33,1	30,0 ± 3,8	26,2 - 33,8
MCHC	g/L	307 ± 48	259 - 355	320 ± 48	272 - 368	328 ± 48,0	280 - 376
	g/dL	30,7 ± 4,8	25,9 - 35,5	32,0 ± 4,80	27,2 - 36,8	32,8 ± 4,8	28,0 - 37,6
RDW-CV	%	13,3 ± 5,0	8,3 - 18,3	13,2 ± 5,0	8,2 - 18,2	13,1 ± 5,00	8,2 - 18,1
RDW-SD	fL	44,0 ± 10,0	34,0 - 54,0	45,6 ± 10,0	35,6 - 55,6	47,2 ± 10,0	37,2 - 57,2
PLT	10 ⁹ /L	67 ± 25	42 - 92	217 ± 45	172 - 262	415 ± 65	350 - 480
MPV	fL	9,1 ± 3,0	6,1 - 12,1	9,5 ± 3,0	6,52 - 12,52	9,4 ± 3,0	6,4 - 12,4
PDW	%	11,6 ± 3,0	8,6 - 14,6	12,1 ± 3,0	9,1 - 15,0	12,6 ± 3,0	9,6 - 15,6
PCT	%	0,061 ± 0,050	0,011 - 0,111	0,207 ± 0,100	0,107 - 0,307	0,389 ± 0,200	0,189 - 0,589



