

ASSAY VALUES AND EXPECTED RANGES

Vizsgálati értékek és várható tartományok

LOT: 1Z0901

05-03-2027

Instruments: DYMIND DH-31, DH-33, DH-36, DH-36VET ,BIOTA Vabio360							
Parameter		CONTROL "Low"		CONTROL "Normal"		CONTROL "High"	
		LOT: 1Z0901 L		LOT: 1Z0901 N		LOT: 1Z0901 H	
		Mean ± Limit	Range Variation	Mean ± Limit	Range Variation	Mean ± Limit	Range Variation
WBC	10 ⁹ /L	2,48 ± 0,50	1,98 - 2,98	7,53 ± 1,00	6,53 - 8,53	16,95 ± 2,50	14,45 - 19,45
LYMPH%	%	41,3 ± 8,3	33,0 - 49,6	30,8 ± 7,7	23,1 - 38,5	25,2 ± 13,8	11,4 - 39,0
GRAN%	%	50,0 ± 17,5	32,5 - 67,5	61,7 ± 9,3	52,4 - 71,0	69,3 ± 9,0	60,3 - 78,3
MID%	%	8,7 ± 7,4	1,3 - 16,1	7,5 ± 6,5	1,0 - 14,0	5,5 ± 4,7	0,8 - 10,2
LYMPH#	10 ⁹ /L	1,02 ± 0,20	0,82 - 1,22	2,32 ± 0,58	1,74 - 2,90	4,27 ± 2,35	1,92 - 6,62
GRAN#	10 ⁹ /L	1,24 ± 0,43	0,81 - 1,67	4,65 ± 0,70	3,95 - 5,35	11,75 ± 1,76	9,99 - 13,51
MID#	10 ⁹ /L	0,22 ± 0,18	0,04 - 0,40	0,56 ± 0,49	0,07 - 1,05	0,93 ± 0,81	0,12 - 1,74
RBC	10 ¹² /L	2,53 ± 0,20	2,33 - 2,73	4,56 ± 0,30	4,26 - 4,86	5,69 ± 0,40	5,29 - 6,09
HGB	g/L	70 ± 6	64 - 76	134 ± 7	127 - 141	170 ± 8	162 - 178
HCT	%	21,5 ± 3,0	18,5 - 24,5	41,8 ± 4,0	37,8 - 45,8	54,0 ± 4,5	49,5 - 58,5
MCV	fL	84,9 ± 6,0	78,9 - 90,9	91,6 ± 6,0	85,6 - 97,6	94,9 ± 6,0	88,9 - 100,9
MCH	pg	27,5 ± 3,8	23,7 - 31,3	29,3 ± 3,8	25,5 - 33,1	29,9 ± 3,8	26,1 - 33,7
MCHC	g/L	324 ± 48	276 - 372	320 ± 48	272 - 368	315 ± 48	267 - 363
RDW-CV	%	17,1 ± 5,0	12,1 - 22,1	17,2 ± 5,0	12,2 - 22,2	17,7 ± 5,0	12,7 - 22,7
RDW-SD	fL	48,6 ± 10,0	38,6 - 58,6	52,5 ± 10,0	42,5 - 62,5	56,0 ± 10,0	46,0 - 66,0
Plt	10 ⁹ /L	74 ± 25	49 - 99	221 ± 45	176 - 266	413 ± 65	348 - 478
MPV	fL	9,3 ± 3,0	6,3 - 12,3	9,6 ± 3,0	6,6 - 12,6	9,7 ± 3,0	6,7 - 12,7
PDW	%	17,1 ± 3,0	14,1 - 20,1	17,3 ± 3,0	14,3 - 20,3	17,4 ± 3,0	14,4 - 20,4
PCT	%	0,068 ± 0,050	0,018 - 0,118	0,213 ± 0,100	0,113 - 0,313	0,401 ± 0,200	0,201 - 0,601

DIAGON Ltd.

48-52., Baross Str.

Budapest, Hungary

H-1047

